

POLY-SPECTRUM-AM

Long-term ECG Ambulatory Monitoring System



2 or 3-channel Holter monitoring solution



wireless connection with clinician's workstation



one button acquisition start



fast data downloading to computer



software package for long-term ECG monitoring with PQ, QT, HRV, ABP, ST, HRT modules




ECG

• PQ • QT • HRV
• ABR • ST • HRT



Neurosoft



3 CHANNELS —
ROUTINE PRACTICE IN
HOLTER MONITORING



LARGE NUMBER OF
RECORDERS ON ONE
WORKSTATION



DUAL-MONITOR
OPERATION MODE



7 DETECTION AND
ANALYSIS ALGORITHMS:
PQ, QT, ST, HRV, ABP, HRT,
AND ARRHYTHMIAS

Poly-Spectrum.NET software performs the automatic analysis of long-term ECG record to detect QRS complexes and rhythm abnormalities dividing them to ventricular and supraventricular ones and pauses, to measure ST segment elevation/depression, and to calculate heart rate (HR) with maximum possible accuracy! We verified our algorithms using the recommended ECG databases (AHA, MIT-BIH, ROHMINE*) and obtained excellent results:

	SENSITIVITY, %	SPECIFICITY, %
ROHMINE	99,98	99,98
MIT-BIH	99,92	99,78
AHA	99,96	99,89

* Russian Society of Holter Monitoring and Non-invasive Electrophysiology



POLY-SPECTRUM-AM



PATIENT BUTTON




BLOOD PRESSURE MONITOR CONNECTION



DEVICE STATE INDICATION



2 BATTERIES OF AA TYPE



CABLES FOR 2- AND 3-CHANNEL ACQUISITION



RESPIRATORY CHANNEL



ACTIVITY AND BODY POSITION CHANNELS




PACEMAKER PULSE DETECTION

CLINICAL APPLICATION


- ▶ heart rate and conduction abnormalities
- ▶ evaluation of antiarrhythmic therapy efficacy
- ▶ study of the autonomic regulation (heart rate variability)
- ▶ detection of ischemic changes
- ▶ evaluation of sudden cardiac death risk
- ▶ assessment of cardiac pacemaker functioning
- ▶ evaluation of cardiac activity in children and adolescents

POLY-SPECTRUM-AM.NET FEATURES



Record Preview

To obtain the higher accuracy of ECG findings during the automatic analysis, you can review the record and exclude the artifacts from it.



Clusters of ECG Complexes


During the analysis the identified QRS complexes are clustered together. Just view two dozens of clusters and get impression on hundred thousand of QRS complexes. Each QRS complex and cluster can be edited manually! The battery of integrated on-screen tools allows both combining, dividing and removing clusters and also arranging the QRS complexes by their similarity. Even the visualization of normal and premature QRS complexes is color-coded that streamlines the assessment.

These advanced features allow clinician to select precisely the most representative clinical ECG data and make a fast and accurate diagnosis.



HRV Analysis

The results of time domain and frequency-domain analyses of heart rate variability (HRV) are presented as spectrograms, bar charts and scattergrams. The time-frequency representation of power spectrum shows the autonomic changes during 24 hours. Most of the graphs and trends are live and you can capture on-the-fly the related QRS complex on the smart ribbon.



ABP Analysis

The software can simultaneously analyze the obtained ECG and ABP data. The ABP results are recorded with BPLab system. The ABP findings can be compared with ECG ones using ABP trends, tables with BP values and loads.



**SUPPORT
OF 12-CHANNEL
RECORDER**




Smart Ribbon to Review and Edit the Record and to Customize the Event Filters

The smart ribbon is a forward-thinking tool to review and edit the QRS complexes and arrhythmia events. You can arrange the markers, define the QRS complex type and morphology, and also review all QRS complexes clustered with user or automatic filters in fast scroll mode. You can customize the event filters and also create your own ones using the multiple parameters.



Analysis of Arrhythmias, PQ, QT, ST, HRT

Using the arrhythmia analysis window you can review the heart rate statistics, including the one for each rhythm abnormality. The data on "PQ analysis", "QT analysis", "ST analysis" and "HRT" pages is clearly shown as trends, graphs and tables. The specially developed tool ensures the color-coded comparison of epochs and also matches them with the reference values. All obtained findings are interactive. Position the mouse cursor to the particular point on graph and you are automatically switched to the related QRS complex.








Automatic Report Generation

There is no need to arrange manually the data obtained during the analysis anymore! Poly-Spectrum-AM.NET software generates the report automatically. The report template can be adjusted according to your needs. It can include the patient information, exam date, record parameters and also tables with obtained data, graphs, ECG fragments, arrhythmias, HR trend, results of ST, HRV, QT and HRT analysis. The automatically generated interpretation can be checked and edited by the clinician with the use of integrated glossary.

ECG PRODUCT LINE



	ECG Leads	Transesophageal Leads	Automatic Interpretation	Detection of Cardiac Pacemaker Pulses	Interface	Application
 Poly-Spectrum-8	12	-	-	+	Wire: USB	Routine ECG test
 Poly-Spectrum-8/E	12	-	+	+	Wire: USB	Routine ECG test with measurement and interpretation, QT, HR, PWV
 Poly-Spectrum-8/EX	12	-	+	-	Wireless: Bluetooth	Routine ECG test with measurement and interpretation, evaluation of autonomic status (HRV), stress ECG, QT, HR, PWV
 Poly-Spectrum-8/G	12	-	+	-	Wireless: Bluetooth, GSM	Telemedicine
 Poly-Spectrum-12/E	14	2	+	+	Wire: USB	Routine ECG test with measurement and interpretation, transesophageal ECG, stress ECG, QT, PWV



Neurosoft

www.neurosoft.com, info@neurosoft.com
 Phones: +7 4932 24-04-34, +7 4932 95-99-99
 Fax: +7 4932 24-04-35
 5, Voronin str., Ivanovo, 153032, Russia

October
2017