
NEURO-AUDIO

Clinical ABR & OAE Analyzer, Screening Audiometer

US-EDITION



NEURO-AUDIO NEW STANDARD IN AEP

Neuro-Audio is the most comprehensive AEP diagnostic device on the market. It is also a diagnostic OAE device and a screening audiometer. It is designed to meet the most stringent requirements of today and is ready for the future. It can be easily adapted for daily clinical routine (including newborn hearing screening and follow-up diagnostics) or for advanced cutting-edge research.



WIDE RANGE OF FEATURES

IN A COMPACT UNIT

Neuro-Audio has a high-quality two-channel amplifier and a built-in auditory stimulator with a wide range of intensities and many supported transducers. It is a portable device that has exceptionally high quality of recording. It can record non-sedated ABR easily and does not require a shielded room or audiometric booth!



OAE DPOAE, TEOAE, SOAE



ASSR
ABR, eABR
MLR, LLR (CAEP), ECochG
Cognitive ERP (P300, MMN)







"Working with the Neuro-Audio system is working with a very user-friendly system. The system is robust and open, which creates many opportunities in order to obtain reliable responses. The service and willingness of the people from Neurosoft to improve their system is certainly an asset".

Robby Vanspauwen, PhD, Clinical Scientist, Dept. for Otorhinolaryngology, Sint-Augustinus Hospital, European Institute for ORL, Antwerp, Belgium

POWER OF INTEGRATION

OPTIMIZED FOR YOUR FEFICIENCY

The common database and user interface, built-in protocols for all tests, integrated reference values and automatic algorithms – all of these are ultimately intended to increase your efficiency and optimize test time usage.

FLEXIBLE SETTINGS

You can adjust almost any software setting and device parameter the way you want it and save it as a new protocol for repeated use.

RELIABLE STORAGE OF PATIENT RECORDS AND TEST DATA

All your patient records and test data are stored safely and securely in a single database used by all other Neurosoft devices. Autosave and backup features protect against data loss.

CUSTOMIZED REPORTS

You can customize design and content of your reports using a powerful report template editor. Also you can use a built-in word processor for advanced editing of your reports.

DATA FXPORT

You can easily export the data and PDF reports to another computer or integrate them with your hospital information system (GDT, HL7 protocols are supported).

OUR SOFTWARE SPEAKS YOUR LANGUAGE

The whole graphical user interface and printed reports are translated to many of the world's languages. This makes our software significantly easier to use no matter which country you are from.

SUPPORTED TRANSDUCERS:

INSERT EARPHONES (ER-3C, ETC.)

AUDIOMETRIC HEADPHONES (TDH-39, ETC.)

BONE VIBRATOR (B-71) LOUDSPEAKERS









ABR

CLINICAL APPLICATION

- Objective audiometry (hearing threshold search using wave V as the indicator).
- Neurology (identification of pathologies in the VIII cranial nerve and auditory pathways of the brainstem).
- Fitting of a cochlear implant (electrically evoked ABR).

CALCULATION OF RESIDUAL NOISE AND FMP (KNOWING WHEN TO STOP)

During recording it's recommended to rely on objective parameters of the response. The residual noise (RN) value helps you determine when to stop averaging. Fmp value indicates presence or absence of a response. Our software calculates these values automatically and shows their reference values. It provides confidence in your results.

WIDE RANGE OF POSSIBILITIES FOR WAVEFORM ANALYSIS

Superimpose ABR waveforms for visual assessment of reproducibility. Just place markers to see the latencies, intervals and amplitudes in analysis tables. Compare results with reference values on the latency/intensity chart. This will improve your efficiency.

OPTIMIZED STIMULUS (CHIRP)

The Chirp and Chirp-LS stimuli are designed with a cochlea model in mind (frequency rises with time) to maximize the evoked response. That's why they are optimal for hearing screening and hearing threshold search (wave V is up to twice the size of a click response). Chirp-LS is optimized to be efficient at any stimulus intensity, so it saves your time.

WEIGHTED AVERAGING FOR HIGH-QUALITY RECORDING


The weighted averaging maintains the steadily high record quality even in the condition of a noisy patient (for example, a lot of movements and muscle activity). As a result the waveform morphology is improved and the residual noise is decreased. It allows saving your time and recording non-sedated ABR from young children.

MINIMIZATION OF MAINS INTERFERENCE

The unique adaptive notch filter (and high harmonic filter) together with the "minimize interference" feature (stimulation rate "jitter") allow recording clear waveforms even with low-quality mains and absence of a shielded room. It saves you time and money.

AUTOMATIC HEARING THRESHOLD SEARCH


Our software has automatic hearing threshold search protocol. Also you can manually select the list of intensities used during the test which can be modified on the fly. It increases ease of use.



ABR latency/intensity chart



Comparison of Click and Chirp ABR waveforms



Dual displays with unlimited waveforms

Estimation of behavioral audiograms with ASSR. The results of ASSR testing can be used to estimate the behavioral pure-tone audiogram in cases when traditional audiometric testing can't be performed (for example, in infants and toddlers).



ACCURATE AND OBJECTIVE RESPONSE DETECTION

F-test response detection method accurately and objectively predicts behavioral audiometric thresholds (based on the analysis of the EEG spectrum). The principal component analysis (PCA) method reduces the impact of muscle-induced noises. The fruitful combination of 2-channel recording, Chirp stimulus and weighted averaging allows achieving the most accurate results.

FULL CONTROL

You can independently control any of the 8 frequencies in multi-ASSR test. Depending on patient's state (awake/asleep) you can change the modulation frequency during testing. Also you can change the maximum test duration. Signal (EEG) monitoring is always visible on the screen. All of these allow you to have full control of the test.

MASKING NOISE

Masking (white noise) allows obtaining accurate audiometric thresholds even in patients suffering from unilateral or conductive hearing loss.

AUTOMATIC HEARING THRESHOLD SEARCH


The hearing threshold search is done automatically with a specified step (in dB) and within selected stimulus intensity range. This makes the test significantly easier to perform and saves your time.

RESIDUAL NOISE CALCULATION

The monitoring of residual noise (RN) and response amplitude (A) for each frequency helps answering the question whether to stop the recording or continue it. It gives you confidence in the obtained results and shortens the test time.



CHIRP STIMULUS







ASSR response probability charts

dB nHL	500 Hz	1000 Hz	2000 Hz	4000 Hz
15			19% 65nV	
10		87% 33nV	50% 45mV	
20	84% 41nV	99% 58nV	99% 28nV	
25		99% 111nV	99% 35nV	95% 18nV
35		100% 25nV	99% B6nV	
30	100% 75nV		100% 89nV	99% 15nV
40	99% 94nV	9156-23nV	99% 20nV	99% 41nV
50	99% 81nV	99% 53nV	99% 51nV	99% 25nV

Table of ASSR trials



- Diagnostics of Meniere's disease/endolymphatic hydrops (ECochG).
- Objective test of sound/speech recognition in the auditory pathway from the ear to the cortex before/after hearing aid fitting (CAEP, cognitive ERP).

FULL SPECTRUM OF AEP: FROM DIAGNOSTIC ABR TO SPECIALIZED AEP

Neuro-Audio can be used to record diagnostic auditory brainstem response (ABR) as well as other diagnostic and research auditory evoked potentials (AEP). The software includes predefined protocols for all AEP types in order to increase your performance. If you perform AEP research studies, you can fully customize all the hardware and software settings. You are free to experiment!

* Ferraro, et. al.- Utility of Area Curve Ratio Electrocochleography in Early Meniere Disease; Arch Otolaryngol Head Neck Surg / Vol. 129, May 2003

ELECTROCOCHLEOGRAPHY: CALCULATION OF SP/AP AREA RATIO

Besides SP/AP amplitude ratio calculation, Neuro-Audio software calculates SP/AP area ratio. Just place three markers on the response waveform. The area ratio calculation makes ECochG significantly more sensitive to diagnose Meniere's disease (according to John Ferraro*).



FOR DIAGNOSTICS AND RESEARCH



Behavioral pure tone audiometry allows obtaining frequency-specific hearing thresholds (audiograms) for air and bone conduction and also in a sound field. It is used for adults and older children (starting from 5 years of age) who can reliably demonstrate a change in behavior when a test sound is heard.



- Complies with IEC 60645-1 (Type 4: screening/monitoring) requirements
- Air and bone conduction testing, testing in a sound field
- Automatic mode (Hughson & Westlake), manual mode (with mouse and keyboard)
- Contralateral masking noise (white noise)
- Silence mode (audiometric booth is not required)
- High-frequency audiometry (up to 16 kHz) with special headphones


CLINICAL APPLICATION

- Objective analysis of cochlear function for patients of all ages
- Newborn hearing screening



- Complies with IEC 60645-6:2009 (Type 1: diagnostic/clinical) requirements
- Full spectrum of OAE tests: from screening to advanced diagnostics
- DPOAE up to 12 kHz (early objective detection of ototoxic and noise-induced hearing loss)
- Completely automatic test (including probe fitting and in-ear calibration of stimulus)





TEOAE analysis

DPOAE analysis

CAUTION: Federal Law restricts this device to sale by or on the order of a practitioner licensed by the law of the State in which he/she practices to use or order the use of the device Neurosoft www.neurosoft.com, info@neurosoft.com Phones: +7 4932 24-04-34, +7 4932 95-99-99 Fax: +7 4932 24-04-35 5, Voronin str., Ivanovo, 153032, Russia