




www.neurosoft.com, com@neurosoft.ru  
Phones: +7 4932 24-04-34, +7 4932 95-99-99  
Fax: +7 4932 24-04-35  
5, Voronin str., Ivanovo, 153032, Russia

# Magnetic Stimulator **Neuro-MS/D**

Fresh Insight  
into TMS



- Research and clinical applications
  - Monophasic waveform
- Two modes of paired pulse stimulation:  
with one coil and two coils
- 7 inch bright color screen to display  
stimulator parameters



- Full USB control of all stimulator parameters allows deepest integration with EMG and EEG recorders
- Wide range of extensions including rapid biphasic system with up 100 Hz frequency and cooling unit
- Trig in/out sockets for precise synchronization with any external amplifier

# Neuro-MS/D

## Magnetic Stimulator



Single pulse mode




Paired pulse mode\*




## Easy control

- Turn to adjust power or IPI
- Gently press to select parameters



+ 6 different types of coils are supported



\* 2 units are required for paired pulse mode.  
Inter-pulse interval can be set in 0...999 ms range.