

# *Comfortable TMS Chair*



*TMS chair specially developed to ensure comfortable seating with head fixing for long TMS sessions*

Fixed seat height – 65 cm

Multiple independent adjustments of backrest and legrest with 2 motors

Backrest adjustment range – 75° (from -7°)

Legrest adjustment – 35°

Remote control for adjustment

Manually adjustable footrest

Footrest adjustment – 29 cm

Big comfortable armrests

Armrest size – L 57 x W 17 cm


Diameter of twin wheels – 7.5 cm (with total lock)

Weight of chair – 63 kg

Weight capacity – 200 kg

Standard power supply – 230 V, 50-60 Hz (special power supplies upon request)

# *Increased comfort level*



[www.neurosoft.com](http://www.neurosoft.com), [info@neurosoft.com](mailto:info@neurosoft.com)  
Phones: +7 4932 24-04-34, +7 4932 95-99-99  
Fax: +7 4932 24-04-35  
5, Voronin str., Ivanovo, 153032, Russia