

CloudEKG Marketing Promotional Brochure

Compact and Mobile EKG System for Everyday Use

Advanced Palmtop Electrocardiography




CloudEKG Device

CloudEKG is a 12-lead battery-operated device intended for recording and transmitting standard electrocardiogram signals for cardiac monitoring. It collects and transfers the signal for processing to a monitoring PC or handheld device wirelessly using Bluetooth technology, displays and prints EKG signals and exam results. .

Main Specifications:

- Acquiring ECG signals
- Bluetooth Technology
- Displaying, processing and storing ECG signals on a PC or handheld device
- Cardiac monitoring and medical diagnosis support for qualified health care providers



12 Channels EKG in the Palm of your Hand



High Quality Report Printout

TeleEMG, LLC - 65 Arlington Road, Woburn, MA 01801, USA
Toll Free: 1-877-TeleEMG - Fax: 1-781-937-3081 - Email: info@TeleEMG.com